

AUG 30 2000

**510(k) SUMMARY FOR THE
ONCOGENE SCIENCE/Bayer Diagnostics
Total Prostate Specific Antigen (tPSA) Microtiter ELISA**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K002121

1. Submitter Information

Prepared 8/22/00

Sheryl Brown-Shimer
Manager, Clinical Trials
Oncogene Science/Bayer Diagnostics
80 Rogers Street
Cambridge, MA 02142
Phone # 617-492-7289
Fax #617-492-8438

2. Device Information

| | |
|-----------------------------|--|
| Trade Name: | Total PSA Microtiter ELISA |
| Common Name: | PSA Immunoassay |
| Classification Name: | Tumor-Associated Antigen Immunological Test Systems Reclassified to Class II effective 9/19/96 |

3. Predicate Device Information

| | |
|----------------------|---|
| Name: | Bayer Immuno 1™ System Prostate Specific Antigen |
| Manufacturer: | Bayer Corporation |
| PMA Number: | P950021 |

4. DEVICE DESCRIPTION

In prostate cancer patients, serum Total Prostate Specific Antigen (tPSA) levels have been shown to increase or decrease with changes in disease burden. Accurate longitudinal measurements of serum Total PSA during the course of disease and therapy can be used as an adjunctive test in the management of prostate cancer patients. Increases in human serum levels of PSA are observed in prostate cancer, benign prostatic hypertrophy or inflammation of the genitourinary tissues. PSA concentrations are not elevated in serum from patients with cancers of the breast, lung, colon, rectum, stomach, pancreas or thyroid. Longitudinal measurements of serum PSA are valuable for monitoring prostate cancer patients since detectable levels of PSA following radical prostatectomy can indicate disease recurrence while low PSA levels indicate disease-free intervals.

In this 510(k) premarket notification, the performance and clinical safety and effectiveness of the Oncogene Science **Total Prostate Specific Antigen Microtiter ELISA** has been established by comparison to a predicate device, the Bayer Immuno 1™ PSA Assay, in accordance with the "Guidance Document For Submission of Tumor Associated Antigen Premarket Notifications, 510(k), to the FDA". Throughout this submission the assay is also referred to as the Total PSA Manual ELISA. Non-clinical studies indicate this assay is a stable, reproducible, highly specific and sensitive assay in which serum components and therapeutic agents do not interfere. Clinical studies indicate substantial equivalence to the predicate device in patient populations. During this study, Oncogene Science Diagnostics, Inc. was purchased by Bayer Corporation.

The Total Prostate Specific Antigen Manual ELISA is a sandwich enzyme immunoassay which utilizes a mouse monoclonal antibody (MM1) for capture and an alkaline phosphatase labeled polyclonal anti-PSA antibody (MP-2) for detection. Both the capture and detector antibodies specifically bind human PSA, are supplied by our manufacturing site, Bayer Corporation, in Elkhart, Indiana and are the same antibodies used in the predicate device for this study – Bayer Immuno 1™ Total PSA assay. The

capture antibody has been immobilized on the interior surface of the microtiter plate wells. To perform the assay, an appropriate volume of serum is incubated in the coated well to allow binding of PSA to the capture antibody, then the conjugated polyclonal antiserum is reacted with bound PSA after a wash step. The amount of detector antibody bound to PSA is measured with Blue PhosTM substrate catalyzed to produce a colored product, allowing quantitation by spectrophotometry. Six prepared PSA standards (0, 0.5, 1, 5, 10 and 25 ng/ml) allow construction of a standard curve for subsequent quantification of Total PSA in serum samples. BioRad PSA Lyphocheck Immunoassay Plus Controls #370 were used in each run for quality control of assay performance.

5. STATEMENT OF INTENDED USE

The Oncogene Science Total Prostate Specific Antigen Microtiter ELISA is an *in vitro* diagnostic assay intended to quantitatively measure Total PSA in human serum. These PSA values should be used in conjunction with information available from clinical and other diagnostic procedures as an aid in the management of prostate cancer patients.

6. SUMMARY OF STUDIES

Nonclinical performance characteristics of the Total Prostate Specific Antigen Manual ELISA were determined at Oncogene Science. These studies included tests for parallelism, linearity, cross-reactivity, interfering compounds, heterophile and human anti-mouse antibody and rheumatoid factor interference, stability and sterility.

The clinical study was retrospective, using banked serum samples obtained by DOCRO, Inc. from institutional specimen banks. The study was separated into 3 parts, with all samples analyzed in the ELISA at DOCRO, Inc. and all samples analyzed in the predicate device at Oncogene Science.

Part I looked at a monitoring cohort with 4-6 serial serum specimens from 60 patients with prostate carcinoma, obtained from two clinical institutions. Both assay results are compared graphically for progressing, responding and stable patients.

Part II was a Clinical Method Comparison statistically analyzing results from both devices for 193 Normal healthy men and for 300 men with prostate diseases. Cumulative frequency distributions and 90 and 95th order statistics are compared in the normal cohort. Regression analyses determine the functional relationship between the test device and the predicate device for various prostate disease groups. Substantive Equivalence among patient groups for the two devices is shown with mean comparisons.

Finally, Part III determined assay variability and analytical sensitivity by analyzing intra-assay, inter-assay and inter-laboratory variability of three quality control materials, three test control materials and six standards in a 20 day Precision study (NCCLS EP5-T2) performed at three laboratory sites.

7. PERFORMANCE DATA – NONCLINICAL STUDIES

7.1 Antibody and Antigen Characterization.

Both the capture antibody (anti-PSA mouse monoclonal) and the detector antibody (anti-PSA polyclonal) are supplied by our manufacturing site, Bayer Diagnostics in Elkhart, Indiana after purification and characterization including electrophoretic mobility, antibody PI and protein content. The Oncogene Science Total PSA ELISA standards are from Scripps Laboratories, San Diego, CA, CAT# P0724. Two lots of this PSA standard antigen were characterized by non-reducing SDS PAGE and show consistent results across the lots and with literature references.

7.2 Parallelism, Linearity, and Spike and Recovery

When serum samples are serially diluted in Sample Diluent (or diluted to 75%, 50%, 25% and 10% in sample diluent), 90-118% recovery is obtained, validating the Sample Diluent as appropriate for diluting and measuring human serum for PSA. Linearity of the assay was indicated by dilution of patient serum samples with a male serum pool at equally spaced proportions of sample concentration and gave percent recoveries from 93-113%. Analyte spiked into 3 patient sera at 3 spike levels gave

average percent recoveries of expected values (versus sample diluent spiked with analyte) of 97%. Therefore serum sample matrix does not affect the ability of the Oncogene Science Total PSA ELISA to accurately measure Total PSA in serum.

7.3 Cross Reactivity and Interference Testing

Three protease inhibitors, Kallikrein, Trypsin and Chymotrypsin, are in the same family as PSA and show high homology with PSA. When spiked into patient sera at 5 final concentrations (higher than normally present in patient sera), there was no effect on the recovery of PSA from these samples, suggesting that none of these proteins cross-react in the assay.

PSA measurements might be performed while patients are taking vitamins, over-the-counter drugs, or undergoing chemotherapy, therefore, these potential exogenous interferences were spiked into a positive control serum pool which was then analyzed for PSA. Potential endogenous interferences found as common serum components were analyzed similarly. All compounds tested allowed close to 100 % recovery of Total PSA, suggesting that none of these materials interfered with measurement of PSA.

The potential non-specific reactivity of HAMA (human anti-mouse antibodies) and of rheumatoid factor in the ELISA were investigated. Testing a number of HAMA and rheumatoid factor positive serum samples indicated there is no significant interference from such samples in the Oncogene Science Total PSA ELISA.

7.4 High-Dose Hook Effect and End-to-End Variation

A 100 µg/ml stock solution of PSA was diluted 6 times and assayed to show that no high-dose hook effect is observed in the ELISA when samples containing very high levels of Total PSA are assayed in the microtiter ELISA.

For end-to-end variation analysis, timing of critical procedural steps was altered from sample to sample within a single run. These experiments showed that greatly extended incubation times could result in erroneously increased or decreased recoveries

of Total PSA (around 15%). However, during normal use following instructions, expected variance in timing of each of the reagent steps is expected to fall well short of these extended times.

7.5 Plate Coating Variability and Device Stability

When a Serum Control sample at 3.3 ng/ml PSA was assayed over an entire plate, a coefficient of variance of 9% was observed indicating consistent sample recovery. Thus uniform antibody coating is indicated throughout the plate. Total PSA standards exhibit long term stability (>18 months). Performance of Total PSA ELISA kits remains robust and consistent when stored for at least 8-9 months post-manufacture at the recommended storage temperature.

7.6 Reagent Bioburden and Preservative Effectiveness Evaluation

Multiple aliquots of the liquid components of the Oncogene Science Total PSA ELISA (all greater than 8 months post-manufacture) were tested at Mycoscience, Inc. No reagent was capable of maintaining or promoting growth of aerobic bacteria, yeast or mold.

7.7 Sample Handling

PSA recovery was analyzed after storage of patient serum samples or BioRad Lyphocheck Immunoassay Plus Controls at 4°C and -20°C for 32 days, as well as after multiple freeze/thaw cycles prior to assay. PSA is stable in human serum for all of these conditions, including after at least 6 cycles of freezing and thawing.

8. PERFORMANCE DATA – CLINICAL STUDIES


Part I – Monitoring. In 58/60 patients, test and predicate device longitudinal graphs are close to or completely superimposable. In one case, a middle sample is very different and in one case, both devices give very high PSA values but the predicate

device plateaus out at a high 500 ng/ml, while the ELISA continues to increase. In all 60 cases, the two devices show the same trend, indicating substantial equivalence of the ELISA to the predicate device as an aid in the management (monitoring) of prostate cancer patients. Further, in most cases there is a direct relationship between changes in Total PSA concentration and clinical course of disease, demonstrating the clinical utility of this assay, when used in conjunction with other clinical indicators, to confirm response to therapy and to signal possible recurrence of malignant disease.

Part II – Clinical Method Comparison.

A) Cutoff - The cumulative distributions of a normal healthy cohort of 200 men ages 40-49 were similar between the predicate and test devices. Analysis of a revised cohort (seven outliers removed) shows that the 95th order statistic for the predicate device, Immuno 1™, is 3.16 ng/ml PSA (95% confidence interval:2.79-3.69), while for the ELISA it is 3.51 ng/ml PSA (95% CI:3.02-4.05). A Demmings Method Comparison of results from the two assays gives a Slope of 1.098 (Identity Line = 1.0). The predicate and test devices give comparable results on the 193 Normal Patients. The Distribution of PSA is shown in the chart below.

Distribution of ELISA PSA (ng/ml) in Normal Healthy Males 40-49



B) Method Comparison – Regression analysis on a prostate disease cohort of 293 men was used to determine the functional relationship between the test and predicate devices. Based on biopsy results, this group consisted of 29% normal, 18 % BPH, 12 % prostatitis, 5 % PIN/Suspicious and 35 % prostate cancer. Mean Total PSA values were similar in each group except the malignant group, which had higher values. There were no significant differences in PSA values between the predicate and test devices in all biopsy groups by Greenhouse-ANOVA. A Demings method comparison of the 189 men with benign prostate disease shows agreement between the values of the predicate and test devices of:

$$\text{Test device} = 1.04 * (\text{predicate device}) - 0.2$$

This indicates a proportional and constant bias in this sub-cohort of 4% of the predicate device value. This increase is within the %CV of each assay, indicating substantial equivalence of these two devices for analysis of men with benign prostate disease.

Demings comparison for the 104 men from the prostate disease cohort with prostate cancer gave agreement of:

$$\text{Test Device} = 1.01 * (\text{predicate device}) - 0.14$$

This indicates no significant proportional and constant bias, and therefore a substantial equivalence of the two devices for men with prostate cancer.

Finally, for the combined study, 486 samples of both normal and prostate disease patients, in the range 0 to 93.3 ng/ml, the relationship between the Total PSA Microtiter ELISA and the predicate device is described by the equation:

$$\text{ELISA PSA} = 1.0094(\text{predicate device}) + 0.0156$$

$$\text{Correlation coefficient (r)} = 0.98$$

C) Substantive Equivalence Among Patient Groups – Equivalence for each of the three cohorts of Normals, Prostate Cancer, and Benign Prostate Diseases was shown by Greenhouse-Geisser ANOVA. Equivalence was defined as: Mean Total PSA ELISA value of the cohort similar to the Mean Immuno 1™ PSA value.

Expected ELISA Values for Total PSA in Normals and in Disease Groups

| | N | % Distribution of tPSA by Disease Category | | | | MEAN (ng/ml) |
|------------------------|-----|--|-------------------|--------------------|---------------|-----------------|
| | | 0-4.0 ng/ml | 4.1 - 10 ng/ml | 10.1 - 30 ng/ml | > 30 ng/ml | |
| Healthy Normals | 193 | 96 | 4 | 0 | 0 | 1.10 |
| Biopsy Results | 293 | | | | | |
| Normal Biopsy | 85 | 25.9 | 57.6 | 16.5 | 0 | 6.82 |
| BPH | 53 | 24.5 | 54.7 | 20.8 | 0 | 6.94 |
| Prostatitis | 36 | 33.3 | 47.2 | 16.7 | 2.8 | 7.77 |
| PIN/Suspicious | 15 | 20 | 73.3 | 6.7 | 0 | 5.96 |
| Prostate Cancer | 104 | 19.2 | 48 | 27 | 5.8 | 10.96 |
| By Gleason | | | | | | |
| 3 | 1 | 100 | | | | |
| 4 | 6 | 16.7 | 83.3 | | | |
| 5 | 11 | 18.2 | 63.6 | 9.1 | 9.1 | |
| 6 | 56 | 25 | 48.2 | 25 | 1.8 | |
| 7 | 18 | 11.1 | 38.9 | 38.9 | 11.1 | |
| 8 | 10 | 0 | 30 | 60 | 10 | |
| 9 | 2 | 0 | 50 | 0 | 50 | |

9. CONCLUSIONS**9.1 Device Performance**

The Oncogene Science Total PSA Microtiter ELISA is reproducible, shows good linearity and parallelism, and has no cross-reactivity, high-dose hook effect, or interference problems.

9.2 Substantial Equivalence

Comparison of clinical sample results from the Oncogene Science Total PSA Microtiter ELISA and from the Bayer Immuno 1™ PSA, for which there is an approved PMA, demonstrate that the two devices are equivalent with respect to method performance, clinical utility and device safety and effectiveness. The two devices give

0000030

Oncogene Science/Bayer Corp
Manual Total PSA ELISA

510(k) #K002121 Additional Information
Safety and Effectiveness

equivalent results in analysis of men with benign prostate disease and in men with prostate cancer. Longitudinal measurements of PSA using both assays for monitoring men with prostate cancer gave close to or completely superimposable results graphically.

- CONFIDENTIAL AND PROPRIETARY -



DEPARTMENT OF HEALTH & HUMAN SERVICES

AUG 3 0 2000

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Ms. Sheryl Brown-Shimer
Manager, Clinical Trials
Oncogene Science
Bayer Corporation
80 Rogers Street
Cambridge, Massachusetts 02142-1168

Re: K002121
Trade Name: Oncogene Science Total PSA Microtiter ELISA
Regulatory Class: II
Product Code: LTJ
Dated: July 10, 2000
Received: July 13, 2000

Dear Ms. Brown-Shimer:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.


A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Page 1 of 1

510 (k) Number (if known): K002121

Device Name: Oncogene Science Total PSA Microtiter ELISA

Indications For Use:

The Oncogene Science Total Prostate Specific Antigen (tPSA) Microtiter ELISA is an *in vitro* diagnostic assay intended to quantitatively measure Total PSA in human serum. These PSA values should be used in conjunction with information available from clinical and other diagnostic procedures as an aid in the management of prostate cancer patients.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER
PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒
(per 21 CFR 801.109)

OR

Over-the-counter Use ☐

(Optional Format 1-2-96)


(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K002121